


Acceptable and Not Acceptable Swabs



ACCEPTABLE	NOT ACCEPTABLE
Ear Culture (C EAR) Eye Culture (C EYE) GC Screen (C GC) Genital Culture (C GEN) Grp B Screen (C GSTR) MRSA Screen (C MRSA) Ortho Culture (C TIS2) Ortho ANA Cult (C TIS2AN) Placenta Culture (C PLAN) Strep Screen by Cult (C STR) Rapid Step (Rapidstrep A Test) VRE Screen (C VRE) Would Culture (C WD) Wound ANA Cult (C WDAN) Yersinia Cult (C YER) Vibrio Cult (C VIB) Stool Culture (C STT)	Wet Mount (WM) Pertussis Culture (PCRPERT) C diff Toxin PCR (CDIFFTOX PCR) Flu A/B & RSV PCR _On-Site (FLU & RSV) Giardia Antigen (GIAAN) O&P (Ova and Parasites, Stool) Varicella Zoster Culture (C VZ) Chlamydia Culture (C CLAM)

Note: Swabs are a **suboptimal** specimen for C AFB and C F (AFB or Fungus Cultures), but they will NOT be rejected.