




Laboratory Sciences
of Arizona



Banner Health®



HORIZON
Laboratory

For any specimen collected using a container/swab outside of this list, please call the NOCO Microbiology laboratory at 970-820-4156 for approval.


Aerobic Swab - eSwab Rapid Group A Strep (throat) by PCR Bacteria/Yeast screens and cultures	
BD NP Swab w/ 3mL UTM Rapid Influenza A&B/RSV/COVID19 by PCR Expanded Respiratory Virus Panel by PCR Bordetella pertussis/parapertussis by PCR	
COPAN Swab Prenatal Group B Strep by Amplification	
BD Affirm Collection Device Vaginitis testing (automated wet prep) Screens for: Candida, Trichomonas, Gardnerella	

**Once printed this document may no longer be current.
Please refer to the online Laboratory Services Directory (JDOS) for current information.**



Gen-Probe Aptima Urine specimen male and female for GC/Chlamydia	
Gen-Probe Aptima Unisex Endocervical Male Urethral GC/Chlamydia (Can also be used for alternative site testing for GC/Chlamydia)	
Para-Pak Ova and Parasite Exam *Do not refrigerate	
Para-Pak Stool Pathogens (GI Panel) by PCR Giardia and Cryptosporidium Antigen	



Sterile Container Clostridium Difficile Lactoferrin Urine (Collection Only – Do not send to lab)	
Urine Collection Device Use either device to transfer urine from collection container to BD tubes	
BD No Additive Yellow Top Tube Any Urinalysis Urine Chemistries	
BD Preservative Grey Top Tube Urine Culture	