

95536.5092 Most Common Unidentified Tests

Copy of version 4.0 (approved and current)

**Last Approval or
Periodic Review Completed** 7/21/2026

**Next Periodic Review
Needed On or Before** 7/21/2027

Effective Date 7/21/2026

Controlled Copy ID 599153

Location Lab Test Directory

Organization LGH Main Lab - Lancaster General Hospital

Author

Nikki Humphreys

Approval and Periodic Review Signatures

Type	Description	Date	Version	Performed By	Notes
Approval	Sendouts Supervisor	7/21/2026	4.0	<i>Nicole Humphreys</i> Nikki Humphreys	

Signatures from prior revisions are not listed.

Version History

Version	Status	Type	Date Added	Date Effective	Date Retired
4.0	Approved and Current	Major revision	7/21/2026	7/21/2026	Indefinite

Linked Documents

- 95536.325 Ordering the Unidentified Test
- 95536.583 Unidentified Test (UNID) Ordering Process Job Aid
- 95536.4545 Unidentified Test (UNID) Specimen Collection and Processing Form

Most Common Unidentified Tests

95536.50924.0
7/21/2026

Test Name	Order	Specimen	Reference lab	Reference lab test number	Sample Processing Requirements/Comments
2,3-Dinor-11 Beta-Prostaglandin F2 Alpha, 24-hour urine	UNID	24HR Urine, no preservative	Quest	10109	24 hour urine, no preservative
2,3-Dinor-11 Beta-Prostaglandin F2 Alpha, Random Urine	UNID	Random urine in sterile cup	Quest	10041	Random urine, refrigerate
3-Hydroxy-3-Methylglutaryl-Coenzyme A Reductase (HMGCR) Antibody (IgG)	UNID	Red	Quest	39044	Spin, aliquot serum, refrigerate
Acylcarnitine, Plasma	UNID	Sodium heparin	Quest	14531	Spin, aliquot plasma, freeze immediately
ACTH Stimulation test					This test involves administering a drug for measurement of the hormone. It is the physician office or floor's responsibility to contact centralized scheduling, the pharmacy for the drug administration, and to provide instructions for the timing and place of the blood draws. This is not collected at lab testing locations.
Alkaline Phosphatase, Bone Specific	UNID	Gold	Quest	29498	Spin, refrigerate tube
Alpha-Gal (Galactose-Alpha-1,3-Galactose) IgE	UNID	Gold	Quest	10554	Spin, refrigerate tube
Alpha-Gal Panel (Includes Galactose-Alpha-1,3 Galactose IgE, Pork IgE, Beef IgE, Lamb IgE)	UNID	Gold	Quest	10555	Spin, refrigerate tube
Alpha-Galactosidase, Serum	UNID	Gold	Mayo Clinic	AGAS	Spin, aliquot serum, freeze
ANALyzeR™ ANA, IFA with Reflex Titer/Pattern, Systemic Autoimmune Panel 1	UNID	Gold x3	Quest	36378	Often ordered by Lancaster Family Allergy. Overnight fasting is preferred. Spin, refrigerated tubes.
Beta-hCG, Quantitative (Tumor Marker)	UNID	Gold	Quest	900726	Spin, refrigerate tube
B-Type Natriuretic Peptide (BNP)	UNID	Lavender	Quest	37386	Spin and aliquot plasma within 1 hour of collection. Freeze plasma immediately. Note: this test is different than NT-proBNP testing that is performed in the Core Lab.
Babesia duncani (WA1) Antibody (IgG), IFA	UNID	Red	Quest	17231	Spin, aliquot serum, refrigerate
BRCA Panel (BRCA1, BRCA2)	UNID	Lavender	Quest	91863	Do not spin, refrigerate
Bullous Pemphigoid BP180 Antibody	UNID	Gold	Quest	16034	Spin, refrigerate tube

Bullous Pemphigoid BP230 Antibody	UNID	Gold	Quest	16136	Spin, refrigerate tube
Chromatin (Nucleosomal) Antibody	UNID	Gold	Quest	34088	Spin, refrigerate tube
Class I&II Donor Specific Antibody (For Penn Transplant Patients)	UNID	Red	HUP	no test code	Do not spin, room temperature
Class I&II HLA Antibody Screen (For Penn Transplant Patients)	UNID	Red	HUP	no test code	Do not spin, room temperature
Class I&II HLA Typing (For Penn Transplant Patients)	UNID	Pink	HUP	no test code	Do not spin, room temperature
Clobazam and Metabolite, Serum/Plasma	UNID	Red	Quest	11409	Spin, aliquot serum, refrigerate
Complement C2	UNID	Red	Labcorp	161414	Spin, aliquot serum, refrigerate
Comprehensive Marfan, Loeys-Dietz, Ehlers-Danlos, and Aortopathy Gene Panel, Varies	UNID	Lavender	Mayo Clinic	CAORG	Do not spin, room temperature
Cyclospora and Isospora Examination	UNID	Stool	Quest	10018	10 grams or 10 mL fresh stool preserved in 1 single Total-Fix® transport vial or 10% formalin. Room temperature.
Dexamethasone Suppression Test					This test involves administering a particlar drug measurement of a particular hormone. It is the physician office or floor's responsibility to contact Centralized scheduling and the pharmacy for the drug administration and to provide instructions for the timing of the blood draws that will occur in Medical Outpatient or as an Inpatient. This is not drawn at a lab testing area.
Drug Monitoring, Methylphenidate Metabolite, Quantitative, Urine	UNID	Urine in sterile cup	Quest	39409	Urine in sterile cup, refrigerate
Drug Screen 9 Panel, Serum or Plasma - Immunoassay Screen with Reflex to Mass Spectrometry Confirmation/Quantitation	UNID	Pink x2	ARUP	0092420	Spin, aliquot plasma, refrigerate within 2 hours of collection
Encephalopathy, Autoimmune/Paraneoplastic Evaluation, Serum	UNID	Red x2	Mayo Clinic	ENS2	Spin, aliquot serum, refrigerate
Eslicarbazepine	UNID	Red	Quest	92139	Spin, aliquot serum, refrigerate

Estrogens, Total, Immunoassay	UNID	Gold	Quest	439	Spin, refrigerate tube
Estrone	UNID	Red	Quest	23244	Spin, aliquot serum, refrigerate
Extended Myositis Panel	UNID	Gold x2	ARUP	3018867	Spin within 2 hours of collection, refrigerate
Ganglioside GQ1b Antibody (IgG)	UNID	Gold	Quest	34144	Spin, refrigerate tube. Fasting is required.
Gastrointestinal Pathogen Panel, PCR, Feces	UNID	Culture and Sensitivity Stool Transport Vial (orange top)	Mayo Clinic	GIP	Stool in a Culture and Sensitivity Stool Transport Vial (orange top), room temperature. Call Sendouts if you need a vial. Specimen stability is only 96 hours from collection.
Hepatitis B Core Antibody (Total) with Reflex to IgM	UNID	Gold	Quest	37676	Spin, refrigerate tube
Histamine, Plasma	UNID	Lavender	Quest	36586	Spin, aliquot plasma, freeze
HLA-A29 Determination	UNID	Lavender x2	Quest	12343	Do not spin, room temperature
HLA-B51 Determination	UNID	Lavender x2	Quest	12344	Do not spin, room temperature
HLA Current Sample (For Penn Transplant Patients)	UNID	Red	HUP	no test code	Do not spin, room temperature
Interleukin 2 Receptor, Soluble, Serum	UNID	Red	ARUP	0051529	Spin, pour off serum and freeze within 30 minutes of collection
Interleukin-6 (IL-6), Serum	UNID	Red	ARUP	0051537	Spin, pour off serum and freeze within 30 minutes of collection
Itraconazole level	UNID	Red	Quest	94092	Spin, aliquot serum, refrigerate
Leptin	UNID	Gold	Quest	90367	Spin, refrigerate tube
Leukotriene E4, 24 hour Urine (LTE4)	UNID	24 hr Urine, no preservative	Quest	11976	24 hour urine, no preservative
Litholink					We do not accept LITHOLINK requests. The patient will need to go to a Labcorp collection site.
LRP4 Autoantibody Test	UNID	Gold	Quest	94744	Spin, refrigerate tube
MyoMarker 3 Plus Profile	UNID	Gold x2	Labcorp	520085	Spin, aliquot serum within one hour of collection, refrigerate
NASH FibroSure® Plus	UNID	Gold x2	Labcorp	550960	Patient should be fasting for at least eight hours. Spin, aliquot serum, freeze
Neurofilament Light Chain, Plasma	UNID	Lavender	Quest	13979	Spin, aliquot plasma, refrigerate

Neurofilament Light Chain, Serum	UNID	Gold	Labcorp	140455	Spin, refrigerate tube
Neuron Specific Enolase (NSE), Serum	UNID	Gold	Quest	34476	Spin, refrigerate tube
Organic Acids, Comprehensive, Quantitative, Urine	UNID	Urine	Quest	35819	Collect urine in sterile cup. Freeze immediately.
Poliovirus (Types 1, 3) Antibodies, Neutralization	UNID	Gold	Quest	94106	Spin, refrigerate tube
Posaconazole	UNID	Red	Quest	94010	Spin, aliquot serum, refrigerate
Primidone and Phenobarbital	UNID	Red	Quest	751	Spin, aliquot serum, refrigerate
Prostaglandin D2 (PG D2), Random Urine	UNID	Urine	Quest	10180	Collect urine in sterile cup. Freeze immediately.
Prostate-specific Antigen (PSA) Total, with Reflex to PSA Free	UNID	Gold	Quest	17569	Spin, aliquot serum, freeze
Prostate-specific Antigen (PSA), Ultrasensitive	UNID	Gold	Labcorp	140731	Spin, refrigerate tube
Prostate-Specific Kallikrein, 4Kscore	UNID	Gold x2	ARUP	2014059	Spin, aliquot serum, freeze
Quest-AD-Detect® Phosphorylated tau217 (p-tau217), Plasma	UNID	Lavender	Quest	13825	Spin, aliquot plasma, refrigerate
Retinol Binding Protein	UNID	Gold	Quest	791	Fasting for at least 12 hours is required. Spin, refrigerate tube.
Schistosoma Antibody (IgG), FMI	UNID	Gold	Quest	34306	Spin, refrigerate tube
Streptococcus pneumoniae Antibody (IgG) (23 Serotypes), MAID	UNID	Gold	Quest	16963	Spin, refrigerate tube
Systemic Sclerosis (Scleroderma) 12 Anitbodies Panel 2	UNID	Gold	Quest	94685	Spin, refrigerate tube
Testosterone, Free	UNID	Red	Quest	18944	Spin, aliquot serum, refrigerate
Testosterone, Free, Bioavailable and Total, MS (Aliases: Bioavailable Testosterone,Testosterone, Free & Weakly Bound)	UNID	Red	Quest	14966	Spin , aliquot serum, refrigerate
Thyroglobulin and Anti-Thryoglobulin to USC Endocrine Laboratory	UNID	Gold	USC Endocrine Lab	Tg + TgAb	Spin, refrigerate tube
Tick-borne Disease, Acute Molecular Panel	UNID	Lavender	Quest	94322	Do not spin, refrigerate
Troponin I, High Sensitivity	UNID	Gold	Quest	34693	Spin, aliquot serum, freeze immediately

Ustekinumab Quantitation with Antibodies, Serum	UNID	Gold	Quest	39986	Spin, refrigerate tube
Vascular Endothelial Growth Factor (VEGF)	UNID	Lavender	Quest	14512	Spin, aliquot plasma, freeze
Vedolizumab Quantitation with Antibodies, Serum	UNID	Red	Quest	36517	Spin , aliquot serum, refrigerate
Vitamin B5 (Pantothenic Acid)	UNID	Red (Protect from light)	Quest	91030	Spin, aliquot serum, protect from light, freeze
Vitamin B6 Profile (Pyridoxal 5-Phosphate and Pyridoxic Acid), Plasma	UNID	Sodium heparin (Protect from light)	Mayo Clinic	B6PRO	Fasting: 12 hours, required Spin, aliquot plasma within two hours of collection, freeze.