

95536.2211 Job Aid: Restricted Tests for Lab Testing Areas

Copy of version 13.0 (approved and current)

**Last Approval or
Periodic Review Completed** 6/10/2025

**Next Periodic Review
Needed On or Before** 6/10/2026

Effective Date 6/10/2025

Controlled Copy ID 573004

Location Lab Test Directory

Organization LGH Main Lab - Lancaster General Hospital

Approval and Periodic Review Signatures

Type	Description	Date	Version	Performed By	Notes
Approval	Lead Technologist	6/10/2025	13.0	<i>Nicole Humphreys</i> Nikki Humphreys	

Signatures from prior revisions are not listed.

Approvals and periodic reviews that occurred before this document was added to Document Control may not be listed.

Version History

Version	Status	Type	Date Added	Date Effective	Date Retired
13.0	Approved and Current	Major revision	6/10/2025	6/10/2025	Indefinite

Job Aid: Restricted Tests for Lab Testing Areas

95536.221113.0
6/10/2025

TITLE: RESTRICTED TESTS FOR LAB TESTING AREAS

Tests Drawn at Downtown Outpatient Pavilion (DOP) Lab Only

Carbon Monoxide (COOX)	Deliver to the Core Lab ASAP
Cold Agglutinin (Z0302)	
Cryoglobulin, Qualitative (CRYG) Cryoglobulin, Quantitative (CRYQ) Cryofibrinogen (CRYOF)	Overnight fasting is required. Draw tube, place the filled tube in warmer immediately following collection, then deliver to Core Lab immediately.
Neutrophil Function Oxidative Burst (CGD)	Call Sendouts. There are special instructions for this test. Draw Monday – Thursday only. Must be received in Sendouts by noon. Do not draw day before holiday.
Aspirin Resistance (ASPRE)	Must have 2 Greiner blue tops hand delivered to the lab
Plavix Resistance (P2Y12)	Must have 2 Greiner blue tops hand delivered to the lab
TEG (TGLB, TGLBL, TPLTM)	Must be hand delivered to the lab

Tests Drawn at DOP, SOP, or WBH Lab Testing Only

Antiphospholipid Syndrome Diagnostic Panel (ASYND)	Requires platelet-poor plasma verified centrifuge.
Antithrombin III (Z0779)	Requires platelet-poor plasma verified centrifuge.
Ascorbic Acid / Vitamin C (ASCOR)	Requires -70°C freezer.
F8 Antigen (F8AG)	Requires platelet-poor plasma verified centrifuge.
Factor Assays	Examples: V, VIII, IX, XI, XII, etc. (Factor V Leiden is an exception may be drawn anytime and anywhere.)
Mixing Studies (MXPTI, MXPTT) Lupus Anticoagulant (LUPUS)	Draw 4 Light blue tubes, STAT delivery to the Core Lab. SOP/WBH Lab Testing may also draw this test – send STAT to WBH Lab for processing.
Plasminogen Activity (Z1219)	Requires platelet-poor plasma verified centrifuge.
Protein C Activity (PRCA)	Requires platelet-poor plasma verified centrifuge.
Protein S Activity (PRSA)	Requires platelet-poor plasma verified centrifuge.
Pyruvic Acid (PYRUV)	Specimen must be processed within 10 minutes of draw per special handling instructions.
Ristocetin Cofactor (RISTO)	Requires platelet-poor plasma verified centrifuge.
Dilute Russel Viper Venom (DRVVT or Russel)	Requires platelet-poor plasma verified centrifuge.
Troponin – when ordered in a set to include 0H, 1H and/or 3 H time	If drawn at WBH or Campus - must Send by STAT courier to LGH.
Urine for Eosinophils (UEOS)	For best sample integrity, send to LGH or WBH Lab within 30 min of collection.
Von Willebrand Panel (VWPAN)	Requires platelet-poor plasma verified centrifuge.

TITLE: RESTRICTED TESTS FOR LAB TESTING AREAS

Von Willebrand Factor Multimer
(ZG476)

Requires platelet-poor plasma verified centrifuge.

Can be drawn at any Lab Testing location but has special requirements

QNatal Advanced (QNATL)	10ml whole blood collected in Streck cell free tubes (black/tan tiger tops). Streck tubes can be obtained from Sendouts department.
GTT (2 or 3 hour) can be drawn at most lab testing sites.	Contact Core lab with any concerns.
FCI – Flow Cytometry	Draw Monday – Thursday all day; Friday before 1000.
Immune Cell Function (Z1132)	Draw Monday-Thursday. Specimen must be received in Sendouts by noon. Testing sent to Quest Diagnostics.
Rituxan Sensitivity CD20 (RITUX)	Draw specimens Monday – Thursday only. Specimen must be received in Sendouts by noon.
Myoglobin, Urine (Z0474)	Urine must be placed into Special Myoglobin Transport tube immediately. Obtain tubes from Sendouts.
T-SPOT TB (TSPTB)	9mL Greiner lithium heparin no gel tube requires for test. Tubes can be obtained from Sendouts department. Monday – Friday collection only.
Molecular Panel (MOLEP) FISH Panel (FISHP) JAK2 JAK2 Exon 12-13 Mutation Analysis (EXO12) High Sensitivity PNH Evaluation (PNH)	Collect Monday – Friday only. Sample must be received in Sendouts by 1400.
CD4 Lymphocyte Subset (CD4LS) TBNK Lymphocyte Subset (TBNKS)	Collect Monday – Thursday only.

Deliver to Main Lab at LGH 3rd Floor

Semen Analysis

Call Centralized Scheduling 717-544-5941