



Origination Date 5/19/2021
Last Approved 2/17/2025
Version Effective Date 2/17/2025
Last Revised 2/8/2023
Next Review 2/17/2027

Owner Patricia Maypole:
Supv Laboratory
Document Area Lab - Phlebotomy
Applicability Piedmont Athens
Regional
Document Laboratory
Groups

Lab Specimen Collection Manual – Order of Draw

POLICY:

Blood specimens are collected using established Order-of-Draw guidelines.

PROCEDURE:

Use correct order-of-draw/fill to help prevent cross contamination of specimens.

For Vacutainer Tubes: Order-of-draw/fill

- Blood Culture
- Red (Plain / No-Additive)
- Blue (Sodium Citrate)
- Gold (SST Serum Separator)
- Green (Lithium Heparin)
- EDTA (Purple / Pink)
- Grey (Sodium Fluoride)
- other additive tubes

For Microtainer tubes (Bullets): Order of Draw/fill

- Purple (EDTA)
- Green (Lithium Heparin)
- Gold (SST Serum Separator)

All Revision Dates

2/8/2023, 9/3/2021, 5/19/2021

Approval Signatures

Step Description	Approver	Date
	Raford Rogers: MD	2/17/2025

Applicability

Piedmont Athens Regional

COPY