



Origination Date	11/29/2007	Owner	Patricia Maypole: Supv Laboratory
Last Approved	2/17/2025	Document Area	Lab - Phlebotomy
Version Effective Date	2/17/2025	Applicability	Piedmont Athens Regional
Last Revised	2/8/2023	Document Groups	Laboratory
Next Review	2/17/2027		

Lab Specimen Collection Manual – Urine 24hr Collection Instructions

Patient Instructions for "24 Hour Urine Collection"

1. At the beginning of the collection, empty your bladder. Discard this first sample into the toilet and flush (do not use this first sample). **Make a note of the time.**
2. From this point on, collect all urine from each voiding for the next 24 hours. It is very important to save ALL urine (from this point on). For each collection, be sure to void into a clean smaller container first, and then carefully pour it over into the large 24 hour urine container. Keep the large 24 hour urine container refrigerated or on ice at all times, during the collection process.
3. At exactly 24 hours from the time the first urine was voided, empty your bladder into a clean smaller container, and then pour this last urine sample into the large 24 hour urine container. This ends the 24 hour collection process.
4. Return to the lab with the completed 24 hour urine collection and your Physician's Orders for this 24 hour test (**preferably between 7am and 5pm**) as soon as possible at the end of the 24 hours.
5. Container should be labeled with minimum of 2 acceptable identifiers. The two identifiers should consist of Patient Name and (at least) one of the following: Patient's Date of Birth, Patient ID/Medical Record #. Place label on Container (not Lid).
6. If using electronic handheld device, scan labeled container to complete collection process in the computer system. If not using electronic handheld device, enter collection in computer system to complete the collection process.

All Revision Dates

2/8/2023, 1/14/2020

Approval Signatures

Step Description

Approver

Date

Raford Rogers: MD

2/17/2025

Applicability

Piedmont Athens Regional

COPY