

95536.1241 Semen Collection Requirements (English)

Copy of version 2.0 (approved and current)

Last Approval or Periodic Review Completed

11/7/2025

Next Periodic Review Needed On or Before

11/7/2026

Effective Date

10/22/2019

Description

Changed to Form for distribution.

Controlled Copy ID 586854

Location

Test Directory

Organization

LGH Main Lab - Lancaster General

Hospital

Approval and Periodic Review Signatures

Туре	Description	Date	Version	Performed By	Notes
Periodic review	Core Lab Manager	11/7/2025	2.0	Lisa Rossi	
Periodic review	Core Lab Manager	11/29/2024	2.0	Lisa Rossi	
Periodic review	Core Lab Manager	1/26/2024	2.0	Lisa Rossi	
Periodic review	Core Lab Manager	4/18/2023	2.0	Lisa Rossi	
Periodic review	Core Lab Manager	5/25/2022	2.0	Lisa Rossi	
Periodic review	Core Lab Manager	7/22/2021	2.0	Lisa Rossi	
Periodic review	Core Lab Manager	10/14/2020	2.0	Lisa Rossi	
Approval	Lab Director	10/19/2019	2.0	Dr. Bruce King	

Bruce King

Approval Core Lab Manager 10/10/2019 2.0

Lisa Rossi

Signatures from prior revisions are not listed.

Approvals and periodic reviews that occurred before this document was added to Document Control may not be listed.

Prior History

CLFF959D.02

Version History

Version	Status	Туре	Date Added	Date Effective	Date Retired
2.0	Approved and Current	Major revision	10/9/2019	10/22/2019	Indefinite

Linked Documents

• 95536.67 Semen Analysis

Semen Collection Requirements (English)

95536.12412.0 10/22/2019

TITLE: SEMEN COLLECTION REQUIREMENTS (ENGLISH)

SPECIMEN REQUIREMENTS:

The patient must adhere to the following collection requirements to provide an appropriate and complete specimen.

- 1. Specimen should be collected following a 2 day period of abstinence, but do not abstain for more than 7 days.
- Specimen MUST be collected DIRECTLY into a clean, dry, leak- proof container. Avoid lubricants. DO NOT collect in a condom / or transfer semen from a condom into another container.
- 3. **ALL** of the semen must be collected. (The first portion of the specimen has the highest sperm concentration).
- 4. Label the container with patient's name, D.O.B. and date/time of collection.
- 5. The following information **MUST** be provided by the patient on an accompanying form:
 - a. Method of collection
 - b. Type of specimen container
 - c. Days of abstinence
 - d. Collection or transport problems
- 6. Transport/delivery to laboratory should be made within 30 minutes of collection (no longer than 60). During transport, protect the specimen from extremes in temperature and light. Recommendation is to put container in a brown bag and keep warm.
- 7. Deliver to third floor specimen receiving window, along with doctor's written request.

ALL SEMEN ANALYSIS TESTS MUST BE SCHEDULED BY CALLING 544-5941.

IF ANYONE DELIVERS A SPECIMEN TO THE LGH LAB WITHOUT AN APPOINTMENT, IT WILL BE REJECTED. THE DOCTORS REQUEST AND THE MANDATORY INFORMATION FROM THE PATIENT MUST ALSO BE WITH THE SPECIMEN.